

Vacuum Sealers

AN5026 Vacuum Sealer

- ▶ **Self-contained**
No shop air or vacuum source needed
- ▶ **Quick and quiet**
Pneumatic cylinders

The Andersen Vacuum Sealer is self-contained, easy to maintain and allows multiple configurations for a wide range of applications.

- ▶ **Safe**
Dual safety switches protect the operator's fingers during sealing

Emergency **STOP** switch on front panel

- ▶ **Versatile**
9 preset heating and cooling settings easily changed by the user

Appropriate for polyethylene and nylon bags of any thickness

- ▶ **Easy to use**
Digital display and digital control
- ▶ **Easy to maintain**
Automatic detection of broken switch
- ▶ **Competitively priced**



Shown:
Model AN5026
with stand

Validation Services Available

- Testing protocol to validate seals
- Burst testing of packages
- Parameters for proper sealing of bag materials
- Selection of appropriate materials


Andersen AN5026 Vacuum Sealer Specifications

Seal length of 15" or 26"

120V and 230V models available

Sealer wires available in 1/8", 1/4" and 3/8" widths

May be connected to a gas flush, as well as to a central vacuum source, if desired




MODEL 5026

Dimensions:
31.5"W x 13.5"D x 12"H
(80.01cm x 34.29cm x 30.5cm)

Power Usage:
120V +/- 10% 50/60Hz 5.7 amps
230V +/- 10% 50/60Hz 3.0 amps

Stand Dimensions:
30.5"W x 28.5"D x 33.5"H
(77.47cm x 72.77cm x 85.9cm)

Extending shelf:
12"W x 16.5"L
(30.48cm x 41.91cm)



ANPRO ANDERSEN PRODUCTS
sterilization redefined

www.anpro.com
Andersen Products, Inc.
Health Science Park • 3202 Caroline Drive
Haw River, NC 27258 USA
800-523-1276 • fax 336-376-8153

www.andcal.com
H. W. Andersen Products of California, Inc.
Health Science Park • 3151 Caroline Drive
Haw River, NC 27258 USA
800-524-3455 • fax 336-376-3088

Distributed in Europe by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea
Essex CO15 4XA UK