

Heat Sealer

Andersen Gas Diffusion Technology



Andersen EOGas[®] System sterilization bags require a heavy duty sealer to ensure good, consistent results.

The Andersen AN1091 heat sealer is a reliable, professional quality sealer that provides a secure seal for EOGas sterilization bags. This rugged sealer features adjustable independent seal and cooling timers, and a built-in cutter knife.

Sealer Specifications

- ▶ Seal length: 24 1/2" Seal width: 1/8"
- ▶ Dimensions: 29"W x 6"H x 7"D
- ▶ 115 Volt standard, 220 Volt optional
- ▶ Power(watts): 575, 16 amps.
- ▶ Shipping weight: 47 lbs.
- ▶ 1 year warranty



*AN1091
shown with
optional stand*

Optional AN5099 Sealer Stand is made of heavy gauge stainless steel, and is designed to support heavy EOGas bags without tipping.

Sealer Stand features a pull out shelf to hold EOGas bags during sealing.

Stand Dimensions:

30.5"W x 28.5"D x 33.5"H
(77.47cm x 77.47cm x 85.9cm)

Extending shelf:

12"W x 16.5"L
(30.48cm x 41.91cm)



sterilization redefined

www.anpro.com
Andersen Products, Inc.
Health Science Park • 3202 Caroline Drive
Haw River, NC 27258-9564 USA
800-523-1276 • fax 336-376-8153

www.andcal.com
H. W. Andersen Products of California, Inc.
Health Science Park • 3151 Caroline Drive
Haw River, NC 27258-9789 USA
800-524-3455 • fax 336-376-3088

Distributed in Europe by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea
Essex CO15 4XA UK